

Biventricular hypertrophic cardiomyopathy in a baby diagnosed with Pompe disease

Pompe hastalığı olan bebekte biventriküler hipertrofik kardiyomiopati


Filiz Ekici

Berna Şaylan Çevik

Mehmet Gündüz#

Departments of
Pediatric Cardiology,
#Pediatric Metabolism,
Ankara Dışkapı Children and
Hematology and Oncology Hospital,
Ankara, Turkey

A four-month-old baby was admitted with respiratory distress and difficulty with feeding. She was born at term with a birth weight of 4000 g and had no symptoms after birth. She was the third child of consanguineous parents. She had hypotonia, a grade 2/6 systolic murmur, tachycardia and mild hepatomegaly. Electrocardiography showed a short PR interval, and biventricular hypertrophy. Echocardiography demonstrated a severe form of biventricular hypertrophy, so that both ventricular cavities were nearly obliterated (Figs. and Video*). The left ventricular mass index (389 g/m^2 , Reference Range (RR); $70\text{-}75 \text{ g/m}^2$) and Z score (+12) were significantly increased. Mild left ventricular outflow tract obstruction was also detected. Serum creatine kinase, creatine kinase-MB, lactate dehydrogenase and alanine aminotransferase levels were elevated: A pro-B-type natriuretic peptide level was markedly increased (4526 pg/ml RR; $0\text{-}100 \text{ pg/ml}$). Dried blood spot analysis revealed no detectable acid alpha glycosidase activity, so the baby was diagnosed with Pompe disease.



Figures– Echocardiography revealed a severe form of biventricular hypertrophic cardiomyopathy. Echocardiographic images obtained at the ventricular-short axis view during the systolic (A) and diastolic (B) phases of the cardiac cycles, parasternal long axis (C) and an apical four chamber (D) view. (E) M-mode echocardiography also showed marked septal hypertrophy. *Supplementary video files associated with this case can be found in the online version. Ao: Aorta; IVS: Interventricular septum; LV: Left ventricle; RV: Right ventricle.

The patient was managed with a beta-blocker and enzyme replacement therapy (Myozyme®, Genzyme Corporation, Cambridge, MA, US, 20 mg/kg infusion every 2 weeks). We observed clinical improvement at short term follow up. Despite multidisciplinary and intensive care management, she unfortunately died due to aspiration pneumonia after two months of therapy.